


Pulse Oximetry



Pulse Oximetry

Philips Reusable SpO₂ Sensors

Product #	Description	Packaging	Use with:*
1 M1191B	Adult Finger Glove Sensor Length = 2m Patient size: >110 lb (>50kg) 8-pin connector	1 sensor	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01, M8105AS, 863276, 863278
M1191BL	AAdult Finger Glove Sensor Length = 3m Patient size: >110 lb (>50kg) 8-pin connector		
M1191T	Adult Finger Glove Sensor Length = 0.45m Patient size: >110 lb (>50kg) 9-pin (D-sub) connector		
2 M1192A	Pediatric/Small-Adult Finger Glove Sensor Length = 1.5m Patient size: 33-110 lb (15-50kg) 8-pin connector		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863051, 863052, 863053, 863054, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; M3921A, M3922A, M3923A, M3924A, M3926A, M3927A, M3928A, M3929A; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01, M8105AS, 863276, 863278
M1192T	Pediatric/Small-Adult Finger Glove Sensor Length = 0.45m Patient size: 33-110 lb (15-50kg) 9-pin (D-sub) connector		
3 M1193A	Reusable Neonatal Foot Wrap Sensor Length = 1.5m Patient size: 2.2 lbs. to 8.8 lbs. (1-4kg) 8-pin connector		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Monitors: SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01, M8105AS, 863276, 863278
M1193T	Reusable Neonatal Wrap Sensor Length = 0.90m Patient size: 2.2 lbs. to 8.8 lbs. (1-4kg) 9-pin (D-sub) connector		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Pulse Oximetry



Pulse Oximetry

Philips Reusable SpO₂ Sensors (continued)

Product #	Description	Packaging	Use with:*
7 M1194A	Pediatric/Adult Ear Clip Sensor Length: 1.5m Patient size: >88 lb (>40kg) 8-pin connector	1 sensor	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01, M8105AS, 863276, 863278
8 M1195A	Reusable Infant Finger Glove Sensor Length: 1.5m Patient size: 8.8 lbs. to 33 lbs. (4-15kg) 8-pin connector		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01, M8105AS, 863276, 863278
9 M1196A	Pediatric/Adult Finger Clip Sensor Length: 3m Patient size: >88 lb (>40kg) 8-pin connector		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01, M8105AS, 860335, 860338, 863276, 863278, 863283
M1196S	Reusable Clip Adult SpO ₂ Sensor Length: 2m Patient size: >40 kg (>88 lbs) 8-pin connector		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01, M8105AS, 863276, 863278

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Pulse Oximetry



10



11

Pulse Oximetry

Philips Reusable SpO₂ Sensors (continued)

Product #	Description	Packaging	Use with:*
10 M1196T	Pediatric/Adult Finger Clip Sensor Length: 90cm Patient size: >88 lb (>40kg) 9-pin (D-sub) connector	1 sensor	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863051, 863052, 863053, 863054, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; M3921A, M3922A, M3923A, M3924A, M3926A, M3927A, M3928A, M3929A; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01, M8105AS, 863276, 863278
11 M1627A	Wristband Velcro straps for M1191A/B SpO ₂ sensors		M1191BL, M1192A, M119B, M1191T, M1192T

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Pulse Oximetry



1



2



3



4

Pulse Oximetry

SpO₂ Single Patient Use sensors

Product #	Description	Packaging	Use with:*
1 M1131A	Adult/Pediatric Finger Clip Single-patient Sensor Length = 45cm Patient size: >44 lb. (>20kg) 9-pin D sub connector	20 sensors	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Telemetry: M2601A transmitter, MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; M3921A, M3922A, M3923A, M3924A, M3926A, M3927A, M3928A, M3929A, Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01, M8105AS, 860335, 860338, 863276, 863278
3 M1132A	Infant Single-patient Wrap Sensor Length = 0.90m. Infant finger Patient size: 6.6 lbs – 22 lbs (3-10kg) 9-pin D sub connector		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; M3921A, M3922A, M3923A, M3924A, M3926A, M3927A, M3928A, M3929A, Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A Telemetry: M8105AT with SpO ₂ Opt. A01, M8105AS, 860335, 860338, 863276, 863278
4 M1133A	Length = 90cm Patient sizes: Adult (finger) >88 lbs (>40kg) Neonate: (foot/hand) <6.6 lbs (<3kg) Infant: (big toe/thumb) 22lbs to 44lbs (10-20kg) 9-pin D sub connector		
M1134A	Adhesive-free Neonatal/Infant/Adult Single-patient Wrap Sensor Length = 90cm Patient sizes: Adult: >88 lb (>40kg) Neonate: < 6.6 lbs (<3kg) Infant: 22-44 lbs (10-20kg) 9-pin D sub connector		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) add the following to the Use with for M1134A: Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863264, 863265, 863266; M3921A, M3922A, M3923A, M3924A, M3926A, M3927A, M3928A, M3929A, Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A Telemetry: M8105AT with SpO ₂ Opt. A01, M8105AS, 863276, 863278

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Pulse Oximetry



1



2



3



4



5

Pulse Oximetry

SpO₂ Adapter Cables for Philips Sensors

Product #	Description	Packaging	Use with:*
1 M1900B	SpO ₂ Adapter Cable (12-pin) Adapts 9-pin sensors for 12-pin sockets Length: 9.8 ft (3m)	1 cable	SpO ₂ Module: Legacy HP/Agilent Monitors: 78352C, 78354C, 78834C, Defibrillators: M1722A, M1722B, M1723A, M1723B, M2475B, SpO ₂ Module: M1020A, Fetal Monitors: M1350B, M1350C
2 M1940A	SpO ₂ Adapter Cable (12-pin) Adapts 8-pin sensors for 12-pin sockets Length: 6.6 ft (2m)		SpO ₂ Modules: M1191A, M1192A, M1193A, M1194A, M1195A, and M1196A SpO ₂ sensors. And Philips instruments: 78352A, 78352C, 78354A, 78354C, 78833A, 78833C, Defibrillators, M1020A, M1722A, M1722B, M1723A, M1723B, M2475B, Monitor Modules, Monitors (General), Monitors (Neonatal/Fetal), M1350B, M1350C
3 M1941A	SpO ₂ Extension Cable (8-pin) Adds additional length between 8-pin sensors and 8-pin sockets Length: 6.6 ft (2m)		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
4 M1943A	SpO ₂ Adapter Cable (8-pin) Adapts 9-pin sensors for 8-pin sockets Length: 3.3 ft (1m)		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) SpO ₂ Modules: M1020B, Multi-Measurement Modules: M3000A, M3001A, M3002A; MP2: M8102A; MP5: M8105A, M8105AT with SpO ₂ Opt. A01; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, M8105AS, 860335, 860338, 863276, 863278
M1943AL	SpO ₂ Adapter Cable (8-pin) Adapts 9-pin sensors for 8-pin sockets Length: 9.8 ft (3m)		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) SpO ₂ Modules: M1020B, Multi-Measurement Modules: M3000A, M3001A, M3002A; MP2: M8102A; MP5: M8105A, M8105AT with SpO ₂ Opt. A01; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, M8105AS, 863276, 863278; HeartStart XL+: 861290
5 M1943NL	OxiMax® SpO ₂ Adapter Cable (8-pin) Adapts 9-pin sensors to 8-pin sockets Length: 9.8 ft (3m) Use only with IntelliVue A02 and A04 sockets		M1020B A02, MP 20 Junior, MP20, MP30, MP40, MP50, MP60, MP70, MP90 A02, M3002A, M8102A, M8105AS, M1020B, M1020B A04, M3001A A04, M3001A, M8105A, M8105AT, M8102A SP4, M8105AS SP4, M8105AT SP4, M8105A SP4

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Pulse Oximetry



1



2



3



4



5



6

Pulse Oximetry

Masimo® LNOP® SpO₂ Adapter Cable

Product #	Description	Packaging	Use with:*
1 989803148221	Masimo SET LNCS adapter cable Adapts LNCS sensors to Philips monitors Length: 10ft (3m)	1 cable	SpO2 Module: M1020A, M1020B; Multi-Measurement Modules: M3000A, M3001A, M3002A; MP2: M8102A; MP5: M8105A with SpO2 Opt. A01; SpO2 SET® Module: M1020B Opt. A03; Multi-Measurement SET® Module: M3001A Opt. A03, M3002A Opt. A03, M8105AS, 863276, 863278
989803182411	rainbowSET 4' cable for rainbow; angle	1 cable	M3001AL, M3002A, M8102A, M8105A, M8105AS - opt A05/SP5
2 989803182421	rainbowSET 4' cable for rainbow; strait	1 cable	863283, 863300, 863301, 863302, 863303, 863304, 863320, 863321, 863322 (check option configuration)
3 989803182431	rainbowSET 10' cable for ALL RB measures	1 cable	
4 989803186471	rainbowSET 12' cable for rainbow; strait	1 cable	M3001AL, M3002A, M8102A, M8105A, M8105AS - opt A05/SP5
5 989803186501	rainbowSET 10' cable for LNCS	1 cable	
6 989803186531	rainbowSET 4' cable for LNCS; straight	1 cable	M3001AL, M3002A, M8102A, M8105A, M8105AS - opt A05 863283, 863300, 863301, 863302, 863303, 863304, 863322 (check option configuration)

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Pulse Oximetry



Pulse Oximetry

Disposable adhesive SpO₂ sensors for FAST or OxiMax (single-patient use)

Product #	Description	Packaging	Use with:*
1 M1901B	Neonatal disposable sensor (single-patient use), sterile Patient size: <7 lbs. or >88 lbs (<3kg or >40kg) 9-pin connector Shelf life: 1 year (not available in N. America or Japan)	24 sensors and tapes	862108, 862231, 862439, 863051, 863052, 863053, 863054, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, 863276, 863277, 863278, 863279, 863283, 865350 (US only), 865351, M1020B, M2703A, M2705A, M2705A, M3001A, M3002A, M8102A, M8105A, M8105AT, M8105AS
2 M1902B	Infant disposable sensor (single-patient use), sterile Patient size: 7-44 lbs (3-20kg) 9-pin connector Shelf life: 1 year (not available in N. America or Japan)		861290, 862108, 862231, 862439, 863051, 863052, 863053, 863054, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, 863276, 863277, 863278, 863279, 863283, 865350 (US only), M1020B, M3001A, M3002A, M8102A, M8105A, M8105AT, M8105AS
3 M1903B	Pediatric disposable sensor (single-patient use), sterile Patient size: 22-110 lbs. (10-50kg) 9-pin connector Shelf life: 1 year (not available in N. America or Japan)		861290, 861493, 861494, 861495, 862108, 862231, 862439, 863051, 863052, 863053, 863054, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, 863276, 863277, 863278, 863279, 863283, 865350 (US only), 865351, M1020B, M1350B, M2703A, M2704A, M2705A, M3001A, M3002A, M3535A, M3536A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, M4735A, M8102A, M8105A, M8105AT, M8105AS
4 M1904B	Adult disposable sensor (single-patient use), sterile Patient size: >66 lbs (>30kg) 9-pin connector Shelf life: 1 year (not available in N. America or Japan)		861290, 861493, 861494, 861495, 862108, 862231, 862439, 862474, 862478, 863051, 863052, 863053, 863054, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, 863276, 863277, 863278, 863279, 863283, 865350 (US only), 864351, M1020B, M1350B, M2703A, M2704A, M2705A, M3001A, M3002A, M3535A, M3536A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, M4735A, M8102A, M8105A, M8105AT, M8105AS

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Pulse Oximetry - Cableless Measurement Supplies



Pulse Oximetry – Cableless Measurement Supplies

Mobile CL SpO₂ Sensors & Accessories

Product #	Description	Packaging	Use with:*
989803165921	Mobile CL 20 single patient SpO ₂ sensors	20 sensors	IntelliVue CL SpO ₂ Pod, 865215
① 989803165941	Mobile CL 20 SpO ₂ sensors + cradles (includes wristbands)	20 sensor kits	
② 989803165951	Mobile CL 20 SpO ₂ cradles(includes wristbands)	20 cradle kits	
③ 989803165961	Mobile CL 50 SpO ₂ wristbands	50 wristbands	
989803168861	Mobile CL SpO ₂ Battery Kit Kit includes replacement battery for Intellivue CL SpO ₂ Pod and required accessories	1 battery kit	

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.