


# Temperature Measurement



# Temperature Measurement

| Product # | Description                                                                                                                          | Packaging         | Use with:*                                                                                                                                                                                                                                                                                         |
|-----------|--------------------------------------------------------------------------------------------------------------------------------------|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 21091A  | Disposable sterile skin surface temperature probe<br>Length: 30" (76cm)<br>Shelf life: at least 12 months                            | 1 box = 20 probes | Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)<br>Temperature Module: M1029A Multi-Measurement Modules: M3000A, M3001A, M3002A SureSigns: 863064, 863065, 863066, 863068, 863087, 863088, M3015A, M3015B MP2 Monitor: M8102A, MP5 Monitor: M8105A            |
| 2 M1837A  | Esophageal/Rectal Temperature Probe, 9 FR<br>Length: 36" (91cm), sterilized, continuous monitoring<br>Shelf life: at least 12 months |                   | Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)<br>Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B |
| 3 21090A  | Esophageal/Rectal Temperature Probe, 12FR<br>Length: 36" (91cm) sterilized, continuous monitoring<br>Shelf life: at least 12 months  |                   | Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)<br>Temperature Module: M1029A; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B                                                                    |

\* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

# Temperature Measurement



# Temperature Measurement

## Disposable Esophageal/Stethoscope Temperature Probes (meets ISO 80601-2-56:2009 Standards)

| Product # | Description                                                                                                                               | Packaging | Use with:*                                                                                                                                                                                                                                                                                              |
|-----------|-------------------------------------------------------------------------------------------------------------------------------------------|-----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ① 21093A  | Esophageal/Stethoscope Temperature Probe, 12 FR<br>Length: 20"(51cm) sterilized, continuous monitoring<br>Shelf life: at least 12 months  | 20 probes | Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)<br>Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088;<br>MP2 Monitor: M8102A, MP5 Monitor: M8105A;<br>MultiMeasurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B |
| ② 21094A  | Esophageal/Stethoscope Temperature Probe, 18FR<br>Length: 20" (51cm) sterilized, continuous monitoring<br>Shelf life: at least 12 months  |           |                                                                                                                                                                                                                                                                                                         |
| 21095A    | Esophageal/Stethoscope Temperature Probe, 24 FR<br>Length: 20" (51cm) sterilized, continuous monitoring<br>Shelf life: at least 12 months |           |                                                                                                                                                                                                                                                                                                         |

\* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

# Temperature Measurement



# Temperature Measurement

## Disposable Foley Catheter Temperature Probes *(meets ISO 80601-2-56:2009 Standards)*

| Product # | Description                                                                                                                      | Packaging | Use with:*                                                                                                                                                                                                                      |
|-----------|----------------------------------------------------------------------------------------------------------------------------------|-----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 M2255A  | Foley Catheter Temperature Probe, 14FR<br>Length: 17" (43cm) sterilized, continuous monitoring<br>Shelf life: at least 12 months | 10 probes | Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)<br>Temperature Module: M1029A; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B |
| 2 21096A  | Foley Catheter Temperature Probe, 16FR<br>Length: 17" (43cm) sterilized, continuous monitoring<br>Shelf life: at least 12 months |           |                                                                                                                                                                                                                                 |
| 3 21097A  | Foley Catheter Temperature Probe, 18FR<br>Length: 17" (43cm) sterilized, continuous monitoring<br>Shelf life: at least 12 months |           |                                                                                                                                                                                                                                 |

\* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

# Temperature Measurement



1



2

# Temperature Measurement

## Reusable Extension Adapter Cables for Disposable Temperature Probes *(meets ISO 80601-2-56:2009 Standards)*

| Product #      | Description                                                                                                                                                           | Packaging | Use with:*                                                                                                                                                                                                                                                                                              |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ① 21082A       | Adapter cable: Connects disposable temperature probes to Philips (or HP or Agilent) monitors<br>Accepts 1/8" mini phone plug. (2-pin connector)<br>Length: 10ft (3m)  | 1 cable   | Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)<br>Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088;<br>MP2 Monitor: M8102A, MP5 Monitor: M8105A;<br>MultiMeasurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B |
| 21082B         | Adapter cable: Connects disposable temperature probes to Philips (or HP or Agilent) monitors<br>Accepts 1/8" mini phone plug. (2-pin connector)<br>Length: 5ft (1.5m) |           |                                                                                                                                                                                                                                                                                                         |
| ② 989803162591 | Adapter cable: Connects disposable temperature probes to Series 400 monitors<br>Length: 10ft (3m)                                                                     |           | Series 400 Temperature Monitors<br>with 1/4" phone plug                                                                                                                                                                                                                                                 |
| ② 989803162601 | Adapter cable: Connects disposable temperature probes to Series 400 monitors<br>Length 5ft (1.5m)                                                                     |           |                                                                                                                                                                                                                                                                                                         |

\* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

# Temperature Measurement




# Temperature Measurement

## Reusable Autoclavable Temperature Probes for Continuous Monitoring (*meets ISO 80601-2-56:2009 Standard*)

| Product #      | Description                                                                                                                                           | Packaging | Use with:*                                                                                                                                                                                                                                                                                                             |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 21075A       | Esophageal/Rectal Temperature Probe (autoclavable), 12FR, 2-prong plug, reusable, continuous monitoring<br>Length: 10ft (3m)                          | 1 probe   | Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)<br>Temperature Module: M1029A, C3: 862474, 862478, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B |
| 2 21076A       | Esophageal/Rectal Temperature Probe (autoclavable), 10FR, 2-prong plug, reusable, continuous monitoring<br>Length: 10ft (3m)                          |           | Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)<br>Temperature Module: M1029A, MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B                                                                                        |
| 3 21078A       | Skin Surface Temperature Probe (autoclavable)<br>0.37" (9.5mm) diameter<br>2-prong plug, reusable tape-on, continuous monitoring<br>Length: 10ft (3m) |           | Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)<br>Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; MultiMeasurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B, 860335, 860338      |
| 4 989803162621 | Esophageal/Rectal Temperature Probe (autoclavable), 12FR, 1/4" phone plug, reusable, continuous monitoring<br>Length: 10ft (3m)                       |           | Series 400 Temperature Monitors with 1/4" phone plug                                                                                                                                                                                                                                                                   |

\* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

# Temperature Measurement



# Temperature Measurement

## Reusable Temperature Probes & Disposable Covers for SureSigns VS1/VS3/VS4 Monitors

| Product #      | Description                                                                                                                                                              | Packaging                  | Use with:*                                                   |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|--------------------------------------------------------------|
| ① 989803162631 | Esophageal/Rectal Temperature Probe (autoclavable), 10FR, 1/4" phone plug, reusable, continuous monitoring<br>Length: 10ft (3m)                                          | 1 probe                    | Series 400 Temperature Monitors with 1/4-inch phone plug.    |
| ② 989803162641 | Skin Surface Temperature Probe (autoclavable)<br>0.37" (9.5mm) diameter<br>1/4in phone plug, reusable tape-on, continuous monitoring<br>Length: 10ft (3m)                |                            |                                                              |
| ③ 989803163281 | Adapter cable required to connect reusable temp probes with 1/4" phone plug to legacy HP/Agilent and Philip 2-pin input connectors<br>Length: 10ft (3m)                  | 1 cable                    | Temperature probes 989803162621, 989803162631, 989803162641  |
| ④ 989803163271 | Adapter cable required to connect reusable temp probes with 2-pin connector to Series 400 temperature monitors with 1/4" phone plug input connector<br>Length: 10ft (3m) |                            | Temperature probes 21075A, 21076A, 21078A                    |
| 989803179611   | Genius 2 Probe Cover (CE marked), disposable, non-sterile, 2 yr min shelf life, size: 17" x 12" x 5.5"                                                                   | 2112 probe covers per case | VS3 (863071 #T02 & 863072 #T02)<br>VS4 (863283 #T02); 866149 |
| 989803179381   | Genius 2 Probe Cover (not-CE marked), disposable, non-sterile, 2 yr min shelf life, size: 17" x 12" x 5.5"                                                               |                            |                                                              |

\* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

# Temperature Measurement



# Temperature Measurement

## Reusable temperature probes & covers for SureSigns VM4/VS2/VS3/VS4 and IntelliVue MP5 monitors

| Product #      | Description                                                      | Packaging    | Use with:*                                                                                                                                                          |
|----------------|------------------------------------------------------------------|--------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 989803143381 | Reusable Oral Temp Probe & Well Kit predictive temp monitoring   | 1 probe      | SureSigns VS3: 863070, 863073, 863074, VM4: 863063, VM6: 863064, 863065, VM8: 863066, IntelliVue MP5: M8105A, M8105AT, M8105AS, VS2: 863081, 863082, 863276, 863278 |
| 2 989803143391 | Reusable Rectal Temp Probe & Well Kit predictive temp monitoring |              | SureSigns VS3: 863070, 863073, 863074, VM4: 863063, VM6: 863064, 863065, VM8: 863066, IntelliVue MP5: M8105A, M8105AT, M8105AS, VS2: 863081, 863082                 |
| M4823A         | Temperature Probe Covers                                         | 1,000 covers | 989803143381, 989803143391 temperature probes, 863276, 863278                                                                                                       |
| 989803180831   | Genius 2 tympanic temperature probe, reusable, non-sterile       | 1 probe      | VS3 (863071 #T02 & 863072 #T02)<br>VS4 (863283 #T02); 866149                                                                                                        |

\* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.